





NR2F6 Antibody

Product Code	CSB-PA016058LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10588
Immunogen	Recombinant Human Nuclear receptor subfamily 2 group F member 6 protein (6-194AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear receptor subfamily 2 group F member 6 (V-erbA-related protein 2) (EAR-2), NR2F6, EAR2 ERBAL2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NR2F6

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA016058LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA016058LA01HU at dilution of 1:100

