



NR3C1 Antibody, Biotin conjugated

Product Code	CSB-PA016059LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04150
Immunogen	Recombinant Human Glucocorticoid receptor protein (1-190AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glucocorticoid receptor (GR) (Nuclear receptor subfamily 3 group C member 1), NR3C1, GRL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NR3C1