



NTSR1 Antibody, HRP conjugated

Product Code	CSB-PA016137LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30989
Immunogen	Recombinant Human Neurotensin receptor type 1 protein (1-67AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neurotensin receptor type 1 (NT-R-1) (NTR1) (High-affinity levocabastine-insensitive neurotensin receptor) (NTRH), NTSR1, NTRR
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NTSR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.