



NTSR2 Antibody, Biotin conjugated

Product Code	CSB-PA016138LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95665
Immunogen	Recombinant Human Neurotensin receptor type 2 protein (359-410AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neurotensin receptor type 2 (NT-R-2) (NTR2) (Levocabastine-sensitive neurotensin receptor), NTSR2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NTSR2