



OPRL1 Antibody, FITC conjugated

Product Code	CSB-PA016360LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41146
Immunogen	Recombinant Human Nociceptin receptor protein (323-370AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nociceptin receptor (Kappa-type 3 opioid receptor) (KOR-3) (Orphanin FQ receptor), OPRL1, OOR ORL1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	OPRL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.