



PLD2 Antibody, Biotin conjugated

Product Code	CSB-PA018145LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14939
Immunogen	Recombinant Human Phospholipase D2 protein (1-165AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phospholipase D2 (PLD 2) (hPLD2) (EC 3.1.4.4) (Choline phosphatase 2) (PLD1C) (Phosphatidylcholine-hydrolyzing phospholipase D2), PLD2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PLD2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.