



PRNP Antibody, Biotin conjugated

Product Code	CSB-PA018739LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04156
Immunogen	Recombinant Human Major prion protein (23-64AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major prion protein (PrP) (ASCR) (PrP27-30) (PrP33-35C) (CD antigen CD230), PRNP, ALTPRP PRIP PRP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PRNP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.