



PRTN3 Antibody, Biotin conjugated

Product Code	CSB-PA018830LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24158
Immunogen	Recombinant Human Myeloblastin protein (1-170AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myeloblastin (EC 3.4.21.76) (AGP7) (C-ANCA antigen) (Leukocyte proteinase 3) (PR-3) (PR3) (Neutrophil proteinase 4) (NP-4) (P29) (Wegener autoantigen), PRTN3, MBN
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PRTN3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.