



# PTGIS Antibody, Biotin conjugated

Product Code	CSB-PA018982LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16647
Immunogen	Recombinant Human Prostacyclin synthase protein (113-199AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prostacyclin synthase (EC 5.3.99.4) (Prostaglandin I2 synthase), PTGIS, CYP8 CYP8A1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	PTGIS