



PTGR2 Antibody, Biotin conjugated

Product Code	CSB-PA018984LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N8N7
Immunogen	Recombinant Human Prostaglandin reductase 2 protein (206-282AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prostaglandin reductase 2 (PRG-2) (EC 1.3.1.48) (15-oxoprostaglandin 13-reductase) (Zinc-binding alcohol dehydrogenase domain-containing protein 1), PTGR2, ZADH1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PTGR2