



PTH1R Antibody, FITC conjugated

Product Code	CSB-PA018988LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q03431
Immunogen	Recombinant Human Parathyroid hormone/parathyroid hormone-related peptide receptor protein (485-593AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Parathyroid hormone/parathyroid hormone-related peptide receptor (PTH/PTHrP type I receptor) (PTH/PTHr receptor) (Parathyroid hormone 1 receptor) (PTH1 receptor), PTH1R, PTHR PTHR1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PTH1R