




PTTG1 Antibody

Product Code	CSB-PA019074LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95997
Immunogen	Recombinant Human Securin protein (2-202AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Securin (Esp1-associated protein) (Pituitary tumor-transforming gene 1 protein) (Tumor-transforming protein 1) (hPTTG), PTTG1, EAP1 PTTG TUTR1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PTTG1

Image



Western Blot

Positive WB detected in: HeLa whole cell lysate, Jurkat whole cell lysate

All lanes: PTTG1 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 23 kDa

Observed band size: 23 kDa