



PTTG1 Antibody, Biotin conjugated

Product Code	CSB-PA019074LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95997
Immunogen	Recombinant Human Securin protein (2-202AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Securin (Esp1-associated protein) (Pituitary tumor-transforming gene 1 protein) (Tumor-transforming protein 1) (hPTTG), PTTG1, EAP1 PTTG TUTR1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PTTG1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.