



PTX3 Antibody, HRP conjugated

Product Code	CSB-PA019079LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26022
Immunogen	Recombinant Human Pentraxin-related protein PTX3 protein (99-164AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pentraxin-related protein PTX3 (Pentaxin-related protein PTX3) (Tumor necrosis factor alpha-induced protein 5) (TNF alpha-induced protein 5) (Tumor necrosis factor-inducible gene 14 protein) (TSG-14), PTX3, TNFAIP5 TSG14
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PTX3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.