





RGS1 Antibody

Product Code	CSB-PA019641LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08116
Immunogen	Recombinant Human Regulator of G-protein signaling 1 protein (14-209AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regulator of G-protein signaling 1 (RGS1) (B-cell activation protein BL34) (Early response protein 1R20), RGS1, 1R20 BL34 IER1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RGS1

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA019641LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA019641LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)