



# RGS1 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA019641LB01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q08116
<b>Immunogen</b>	Recombinant Human Regulator of G-protein signaling 1 protein (14-209AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Regulator of G-protein signaling 1 (RGS1) (B-cell activation protein BL34) (Early response protein 1R20), RGS1, 1R20 BL34 IER1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	RGS1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.