



RGS1 Antibody, FITC conjugated

Product Code	CSB-PA019641LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08116
Immunogen	Recombinant Human Regulator of G-protein signaling 1 protein (14-209AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regulator of G-protein signaling 1 (RGS1) (B-cell activation protein BL34) (Early response protein 1R20), RGS1, 1R20 BL34 IER1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RGS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.