



RGS1 Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA019641LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q08116 |
| Immunogen | Recombinant Human Regulator of G-protein signaling 1 protein (14-209AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Regulator of G-protein signaling 1 (RGS1) (B-cell activation protein BL34) (Early response protein 1R20), RGS1, 1R20 BL34 IER1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | RGS1 |