



S1PR2 Antibody, HRP conjugated

Product Code	CSB-PA020651LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95136
Immunogen	Recombinant Human Sphingosine 1-phosphate receptor 2 protein (293-353AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sphingosine 1-phosphate receptor 2 (S1P receptor 2) (S1P2) (Endothelial differentiation G-protein coupled receptor 5) (Sphingosine 1-phosphate receptor Edg-5) (S1P receptor Edg-5), S1PR2, EDG5
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	S1PR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.