




SEPT3 Antibody

Product Code	CSB-PA021026LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UH03
Immunogen	Recombinant Human Neuronal-specific septin-3 protein (238-358AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuronal-specific septin-3, SEPT3, SEP3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SEPT3

Image



Western Blot

Positive WB detected in: Mouse brain tissue


All lanes: SEPT3 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 41, 10 kDa

Observed band size: 41 kDa



Immunohistochemistry of paraffin-embedded human melanoma using CSB-PA021026LA01HU at dilution of 1:100