



SIPA1 Antibody, FITC conjugated

Product Code	CSB-PA021330LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96FS4
Immunogen	Recombinant Human Signal-induced proliferation-associated protein 1 protein (759-927AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Signal-induced proliferation-associated protein 1 (Sipa-1) (GTPase-activating protein Spa-1) (p130 SPA-1), SIPA1, SPA1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SIPA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.