



SLC11A2 Antibody, FITC conjugated

Product Code	CSB-PA021381LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49281
Immunogen	Recombinant Human Natural resistance-associated macrophage protein 2 protein (1-69AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Natural resistance-associated macrophage protein 2 (NRAMP 2) (Divalent cation transporter 1) (Divalent metal transporter 1) (DMT-1) (Solute carrier family 11 member 2), SLC11A2, DCT1 DMT1 NRAMP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SLC11A2