



SLC11A2 Antibody, Biotin conjugated

Product Code	CSB-PA021381LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49281
Immunogen	Recombinant Human Natural resistance-associated macrophage protein 2 protein (1-69AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Natural resistance-associated macrophage protein 2 (NRAMP 2) (Divalent cation transporter 1) (Divalent metal transporter 1) (DMT-1) (Solute carrier family 11 member 2), SLC11A2, DCT1 DMT1 NRAMP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SLC11A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.