



SLC1A2 Antibody, FITC conjugated

Product Code	CSB-PA021433LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43004
Immunogen	Recombinant Human Excitatory amino acid transporter 2 protein (460-574AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Excitatory amino acid transporter 2 (Glutamate/aspartate transporter II) (Sodium-dependent glutamate/aspartate transporter 2) (Solute carrier family 1 member 2), SLC1A2, EAAT2 GLT1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC1A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.