




SLC25A12 Antibody

Product Code	CSB-PA021478LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75746
Immunogen	Recombinant Human Calcium-binding mitochondrial carrier protein Aralar1 protein (202-317AA)
Raised In	Rabbit
Species Reactivity	Human, Rat, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium-binding mitochondrial carrier protein Aralar1 (Mitochondrial aspartate glutamate carrier 1) (Solute carrier family 25 member 12), SLC25A12, ARALAR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC25A12

Image



Western Blot

Positive WB detected in: HEK293 whole cell lysate, Rat heart tissue, Rat kidney tissue, Mouse kidney tissue, Mouse heart tissue


All lanes: SLC25A12 antibody at 3µg/ml

Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 75, 63 kDa

Observed band size: 75 kDa



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA021478LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA021478LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.