



SLC25A12 Antibody, FITC conjugated

Product Code	CSB-PA021478LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75746
Immunogen	Recombinant Human Calcium-binding mitochondrial carrier protein Aralar1 protein (202-317AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium-binding mitochondrial carrier protein Aralar1 (Mitochondrial aspartate glutamate carrier 1) (Solute carrier family 25 member 12), SLC25A12, ARALAR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC25A12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.