



SLC25A12 Antibody, Biotin conjugated

Product Code	CSB-PA021478LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75746
Immunogen	Recombinant Human Calcium-binding mitochondrial carrier protein Aralar1 protein (202-317AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium-binding mitochondrial carrier protein Aralar1 (Mitochondrial aspartate glutamate carrier 1) (Solute carrier family 25 member 12), SLC25A12, ARALAR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC25A12