



STXBP1 Antibody, FITC conjugated

Product Code	CSB-PA022905LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61764
Immunogen	Recombinant Human Syntaxin-binding protein 1 protein (245-594AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Syntaxin-binding protein 1 (MUNC18-1) (N-Sec1) (Protein unc-18 homolog 1) (Unc18-1) (Protein unc-18 homolog A) (Unc-18A) (p67), STXBP1, UNC18A
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	STXBP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.