



TREM2 Antibody, FITC conjugated

Product Code	CSB-PA024405LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZC2
Immunogen	Recombinant Human Triggering receptor expressed on myeloid cells 2 protein (19-174AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Triggering receptor expressed on myeloid cells 2 (TREM-2) (Triggering receptor expressed on monocytes 2), TREM2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TREM2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.