



TREM2 Antibody, Biotin conjugated

Product Code	CSB-PA024405LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZC2
Immunogen	Recombinant Human Triggering receptor expressed on myeloid cells 2 protein (19-174AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Triggering receptor expressed on myeloid cells 2 (TREM-2) (Triggering receptor expressed on monocytes 2), TREM2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TREM2