




TRIM3 Antibody

Product Code	CSB-PA024467LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75382
Immunogen	Recombinant Human Tripartite motif-containing protein 3 protein (285-487AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tripartite motif-containing protein 3 (Brain-expressed RING finger protein) (RING finger protein 22) (RING finger protein 97), TRIM3, BERP RNF22 RNF97
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TRIM3

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA024467LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.