



TRIM3 Antibody, Biotin conjugated

Product Code	CSB-PA024467LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75382
Immunogen	Recombinant Human Tripartite motif-containing protein 3 protein (285-487AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tripartite motif-containing protein 3 (Brain-expressed RING finger protein) (RING finger protein 22) (RING finger protein 97), TRIM3, BERP RNF22 RNF97
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TRIM3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.