



TRIM3 Antibody, Biotin conjugated

Product Code	CSB-PA024467LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75382
Immunogen	Recombinant Human Tripartite motif-containing protein 3 protein (285-487AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tripartite motif-containing protein 3 (Brain-expressed RING finger protein) (RING finger protein 22) (RING finger protein 97), TRIM3, BERP RNF22 RNF97
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TRIM3