




TWIST1 Antibody

Product Code	CSB-PA025358LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15672
Immunogen	Recombinant Human Twist-related protein 1 protein (1-80AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Twist-related protein 1 (Class A basic helix-loop-helix protein 38) (bHLHa38) (H-twist), TWIST1, BHLHA38 TWIST
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TWIST1

Image



Western Blot
 Positive WB detected in: MCF-7 whole cell lysate
 All lanes: TWIST1 antibody at 3µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 21 kDa
 Observed band size: 50 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.