



TWIST1 Antibody, Biotin conjugated

Product Code	CSB-PA025358LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15672
Immunogen	Recombinant Human Twist-related protein 1 protein (1-80AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Twist-related protein 1 (Class A basic helix-loop-helix protein 38) (bHLHa38) (H-twist), TWIST1, BHLHA38 TWIST
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TWIST1