



UBASH3A Antibody, FITC conjugated

Product Code	CSB-PA025427LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57075
Immunogen	Recombinant Human Ubiquitin-associated and SH3 domain-containing protein A protein (287-455AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-associated and SH3 domain-containing protein A (Cbl-interacting protein 4) (CLIP4) (Suppressor of T-cell receptor signaling 2) (STS-2) (T-cell ubiquitin ligand 1) (TULA-1), UBASH3A, STS2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	UBASH3A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.