



VCP Antibody, Biotin conjugated

Product Code	CSB-PA025813LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55072
Immunogen	Recombinant Human Transitional endoplasmic reticulum ATPase protein (704-806AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transitional endoplasmic reticulum ATPase (TER ATPase) (EC 3.6.4.6) (15S Mg(2+)-ATPase p97 subunit) (Valosin-containing protein) (VCP), VCP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	VCP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.