



TNFRSF9 Antibody, HRP conjugated

Product Code	CSB-PA08109B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07011
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 9 protein (24-186AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor receptor superfamily member 9 (4-1BB ligand receptor) (CDw137) (T-cell antigen 4-1BB homolog) (T-cell antigen ILA) (CD antigen CD137), TNFRSF9, CD137 ILA
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TNFRSF9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.