



TNFRSF9 Antibody, HRP conjugated

Product Code	CSB-PA08109B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07011
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 9 protein (24-186AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor receptor superfamily member 9 (4-1BB ligand receptor) (CDw137) (T-cell antigen 4-1BB homolog) (T-cell antigen ILA) (CD antigen CD137), TNFRSF9, CD137 ILA
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TNFRSF9