



PTN Antibody, HRP conjugated

Product Code	CSB-PA10239B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21246
Immunogen	Recombinant Human Pleiotrophin protein (33-168AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pleiotrophin (PTN) (Heparin-binding brain mitogen) (HBBM) (Heparin-binding growth factor 8) (HBGF-8) (Heparin-binding growth-associated molecule) (HB-GAM) (Heparin-binding neurite outgrowth-promoting factor 1) (HBNF-1) (Osteoblast-specific factor 1) (OSF-1), PTN, HBNF1 NEGF1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PTN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.