



LIN28A Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA13619B0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9H9Z2 |
| Immunogen | Recombinant Human Protein lin-28 homolog A protein (1-209AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Protein lin-28 homolog A (Lin-28A) (Zinc finger CCHC domain-containing protein 1), LIN28A, CSDD1 LIN28 ZCCHC1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | LIN28A |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |