



HMMR Antibody, HRP conjugated

Product Code	CSB-PA14189B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75330
Immunogen	Recombinant Human Hyaluronan mediated motility receptor protein (478-653AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hyaluronan mediated motility receptor (Intracellular hyaluronic acid-binding protein) (Receptor for hyaluronan-mediated motility) (CD antigen CD168), HMMR, IHABP RHAMM
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	HMMR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.