



HMMR Antibody, HRP conjugated

Product Code	CSB-PA14189B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75330
Immunogen	Recombinant Human Hyaluronan mediated motility receptor protein (478-653AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hyaluronan mediated motility receptor (Intracellular hyaluronic acid-binding protein) (Receptor for hyaluronan-mediated motility) (CD antigen CD168), HMMR, IHABP RHAMM
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	HMMR