



PRKCB Antibody, HRP conjugated

Product Code	CSB-PA15849B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05771
Immunogen	Recombinant Human Protein kinase C beta type protein (517-643AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein kinase C beta type (PKC-B) (PKC-beta) (EC 2.7.11.13), PRKCB, PKCB PRKCB1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PRKCB