




ADGRB2 Antibody

Product Code	CSB-PA523985LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60241
Immunogen	Recombinant Human Adhesion G protein-coupled receptor B2 protein (1371-1545AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:500-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adhesion G protein-coupled receptor B2 (Brain-specific angiogenesis inhibitor 2), ADGRB2, BAI2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	ADGRB2

Image



Western Blot

Positive WB detected in: Rat brain tissue


All lanes: ADGRB2 antibody at 3.2μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 173, 172, 170 kDa

Observed band size: 173 kDa



Immunofluorescent analysis of HepG2 cells using CSB-PA523985LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)