



# ADGRB2 Antibody, HRP conjugated

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA523985LB01HU   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                  |
| <b>Uniprot No.</b>         | O60241   |
| <b>Immunogen</b>           | Recombinant Human Adhesion G protein-coupled receptor B2 protein (1371-1545AA)                 |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human  |
| <b>Tested Applications</b> | ELISA  |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | HRP  |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4               |
| <b>Purification Method</b> | >95%, Protein G purified   |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Alias</b>               | Adhesion G protein-coupled receptor B2 (Brain-specific angiogenesis inhibitor 2), ADGRB2, BAI2 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Research Area</b>       | Cardiovascular   |
| <b>Target Names</b>        | ADGRB2   |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                    |