



# ADGRB2 Antibody, HRP conjugated

Product Code	CSB-PA523985LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60241
Immunogen	Recombinant Human Adhesion G protein-coupled receptor B2 protein (1371-1545AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adhesion G protein-coupled receptor B2 (Brain-specific angiogenesis inhibitor 2), ADGRB2, BAI2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	ADGRB2