



ECI2 Antibody, HRP conjugated

Product Code	CSB-PA527983LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75521
Immunogen	Recombinant Human Enoyl-CoA delta isomerase 2, mitochondrial protein (119-219AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Enoyl-CoA delta isomerase 2, mitochondrial (EC 5.3.3.8) (DRS-1) (Delta(3),delta(2)-enoyl-CoA isomerase) (D3,D2-enoyl-CoA isomerase) (Diazepam-binding inhibitor-related protein 1) (DBI-related protein 1) (Dodecenoyl-CoA isomerase) (Hepatocellular carcinoma-associated antigen 88) (Peroxisomal 3,2-trans-enoyl-CoA isomerase) (pECI) (Renal carcinoma antigen NY-REN-1), ECI2, DRS1 HCA88 PECl
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	ECI2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.