





ITPR1 Antibody

Product Code	CSB-PA614538LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14643
Immunogen	Recombinant Human Inositol 1,4,5-trisphosphate receptor type 1 protein (1707-1966AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inositol 1,4,5-trisphosphate receptor type 1 (IP3 receptor isoform 1) (IP3R 1) (InsP3R1) (Type 1 inositol 1,4,5-trisphosphate receptor) (Type 1 InsP3 receptor), ITPR1, INSP3R1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ITPR1


Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA614538LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA614538LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA614538LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.