



# ITPR1 Antibody, HRP conjugated

Product Code	CSB-PA614538LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14643
Immunogen	Recombinant Human Inositol 1,4,5-trisphosphate receptor type 1 protein (1707-1966AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inositol 1,4,5-trisphosphate receptor type 1 (IP3 receptor isoform 1) (IP3R 1) (InsP3R1) (Type 1 inositol 1,4,5-trisphosphate receptor) (Type 1 InsP3 receptor), ITPR1, INSP3R1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ITPR1